



## hCG Rapid Laboratory Pregnancy Test

Vizient Contract # LB0812



- Easy to use, minimal training needed.
- One-step, no additional reagents required.
- Proven accuracy of >99%.
- Compact design maximizes storage space in laboratory and Point Of Care locations.
- Includes hCG variance H-hCG and hCGBcf, for the detection of early and mid-stage pregnancy.
- Combo Test, cut off 10 mIU/hCG serum, 20 mIU/hCG urine, can detect  $1 \times 10^6$  pmol/ml level of hCGBcf and/or H-hCG.

